



## 印製耳型10大程序 Ten Procedures for Taking Ear Impressions

- 


1. 耳窺鏡檢查  
Otosopic Examination
- 

2. 將棉球放入耳內保護耳朵  
Insert otoblock into the ears
- 

3. 預備高質素矽膠印模材料  
Prepare the high quality silicone impression material
- 

4. 將印耳型材料混合  
Mix the ear impression material
- 

5. 印耳型時將口部張開  
Open mouth widely when impression material is being injected into the ear
- 

6. 口部保持張開約五分鐘  
Do not close mouth until the impression material has completely cured
- 

7. 將耳型輕輕取出  
Take out the impression gently
- 

8. 檢查耳型是否完整無缺  
Check the impression to make sure it is complete
- 

9. 耳窺鏡檢查，確保耳道清潔  
Do otosopic examination again to ensure ear canal free of impression material
- 

10. 將耳型存放在盒內  
Store the impression into a box